

2010 | MES | SPERM QUALITY ANALYZERS FOR VETERINARY



SQA-V Series

Sperm Quality Analyzers



SQA-V is a high performance analytical device that combines the most up-to-date electro-optic technology, computer algorithms and video microscopy. The test results are displayed on the operational screen and can be printed out and saved in the internal archive.

- Automatic results in 75 seconds
 - Minimal operator experience or training required
 - Automatic dosimetry calculations performed for semen extension
 - On-screen visualization of the semen sample on the video screen of the SQA-V or on a PC display (with counting grid and image freezing) using a standard laboratory slide or an SQA-V capillary
 - Variable optical magnification from x300 to x500
 - A complete semen analysis report can be automatically printed-out
 - QC: self-testing, self-calibrating
 - Reusable, washable testing capillary
 - PC-compatible with special sperm software: import patient data, capture video, save data, full analysis report
-
- **SQA-Vb** for bulls
 - **SQA-Vp** for pigs
 - **SQA-Ve** for equines
 - **SQA-Vt** for turkeys

QwikCheck Gold Series

Sperm Quality Analyzers



QwikCheck Gold is designed for smaller operations that still demand the highest levels of accuracy and precision during semen evaluation. Internal temperature control and an external pre-heating dock ensure consistent readings throughout the testing cycle.

- Fast, simple and accurate semen testing in less than one minute
 - Temperature control
 - Prints test results on a label maker
 - 250 record archive for test results and later printing
 - QC: self-testing, self-calibrating
 - Reusable, washable testing capillary
-
- **QwikCheck Gold Bull** for bulls
 - **QwikCheck Gold Pig** for pigs
 - **QwikCheck Gold Equine** for equines

QwikCheck QC Series

Sperm Quality Analyzers



QwikCheck QC is designed for Vets who inseminate sows and want to QC semen samples. Internal heating insures accuracy and results can be printed to an optional label printer.

- Automatic assessment of semen parameters
 - Fast, simple and accurate semen testing in less than one minute
 - Temperature control
 - Prints test results on a label maker
 - QC: Self-testing, self-calibrating
 - Reusable, washable testing capillary
-
- **QwikCheck QC Bull** for bulls
 - **QwikCheck QC Pig** for pigs
 - **QwikCheck QC Equine** for equines