

MX 950

Metallurgical microscope

- Trinocular microscope with infinite optics
- Quadruple objective nosepiece
- Coaxial fine and coarse adjustment
- LED illumination adjustable 3 V
- Objectives plan achromat: 5x, 10x, 20x, 50x
- Ingenious stand for convenient operation
- Photo- and video documentation



Specifications

MX 950 metallurgical microscope

Standard set

Head	seidentopf type trinocular head
Eyepiece	WF 10x/18
Objective	Plan Achromat: 5x, 10x, 20x, 50x
Nosepiece	quadruple nosepiece
Stage	— double layer mechanical stage 150x140 mm — movement range 75x50 mm
Focusing	coaxial Fine & Coarse Adjustment
Illumination	Kohler illumination, Epi-illuminator with iris aperture diaphragm and iris field diaphragm, 3 V LED illumination, brightness adjustable
Filters	blue, green, yellow, ground
Additionally (on request)	binocular head, eyepieces, micrometer, C-Mount adapter, filters, objectives, digital cameras, software for management of digital albums

Ordering Information

Description

Code

MX 950 metallurgical microscope, standard set

09.0950.03



We reserve the right to change specification without notice.