

Vision Cyto[®] Pap

Automation of cervical cytology



Automated cervical cytology imaging analyzer

Workflow automation and standardization. Improved quality and increased productivity. Time efficiency.

Automatic scanning


Extended and range focus

Creating virtual sample

Automatic pre-classification
by Bethesda system




Works with LBC slides





Reliability and objectivity, confirmed by documentation

Efficient and professional solution for cytological analysis



Automatic pre-classification by Bethesda system




NILM

- Reactive changes


Squamous / Glandular

- ASC-US / AGC
- ASC-H / AGC-FN
- LSIL
- HSIL / AIS
- SCC / ADC


Main features




Automatic scanning and creating virtual sample



Automatic pre-classification by Bethesda system



Quick verification of results (validation required by a medical technologist)



Hints from the pre-set cell atlas



Database for archive management




Education and professional development



Telemedicine and remote consultations




Remote access



Bi-directional communication with LIS


Automated systems for Pap smear analysis

Vision Cyto® Pap Assist



Successive slide loading by operator

Vision Cyto® Pap 4Pro



Multiple slide loading by operator

Vision Cyto® Pap 8Pro



Slide tray loading

Vision Cyto® Pap Ultimate



Loading of up to 200 slides simultaneously

Specification

Main characteristics

	Vision Cyto® Pap Assist	Vision Cyto® Pap 4Pro
Slide capacity	1 slide	up to 4 slides
Slide handling	Manual	Manual
Slide identification	Barcode-labeled: manual	Barcode-labeled: manual
Content	Scanning microscope, PC and monitor, Vision Cyto® Pap software	Scanning microscope, PC and monitor, Vision Cyto® Pap software
	Vision Cyto® Pap 8Pro	Vision Cyto® Pap Ultimate
Slide capacity	up to 8 slides	up to 200 slides
Slide handling	Slide tray	Automated
Slide identification	Barcode-labeled: manual/automated	Barcode-labeled: automated
Content	Scanning microscope, slide tray (2 pcs.), PC and monitor, Vision Cyto® Pap software, barcode reader (option)	Scanning microscope, PC and monitor, touch-screen monitor for control, Vision Cyto® Pap software, Vision Cyto® Pap Remote software, automatic slide loader, built-in barcode reader, barcode printer for slides

Software

Vision Cyto® Pap	<ul style="list-style-type: none">— Scanning slide and creating virtual sample— Automatic pre-classification by Bethesda system— Database of patients, digital samples and analysis results— A pre-set cell atlas— “Shortcut” keys, masks, quick preview, color marks and comments for captured cells in the sample. Drag&Drop feature— Customizable reference guide to generate reports, following your personal requirements— Remote access and network capabilities
------------------	--



* Product images are shown for reference only and final product may differ

Rev 7.2/06.2017 EN



VISION

www.wm-vision.com

Franz-Siegel-Gasse 1
2380 Perchtoldsdorf, Austria
tel.: +43 (1) 804 81 84
fax: +43 (1) 804 81 85
vienna@westmedica.com



www.westmedica.com