

# SlideVision (B) — Digital Approach and Microscopy Automation



- **Efficiency**

Your style of work is changed completely — from manual to automated


- **Simplicity and convenience**

You work with one system only and do not waste time on manual operations with a microscope


- **Automation**

Automated preparation of digital specimen — several visual fields are merged into one


# Specimen's preparation procedure



# Automatic merge of visual fields into one specimen



# Digital approach and microscopy automation



- **Efficiency**

Your style of work is changed completely — from manual to automated

- **Simplicity and convenience**

You work with one system only and do not waste time on manual operations with a microscope

- **Automation**

Automated preparation of digital specimen – several visual fields are merged into one

- **Dialogue navigation**

If you are interested in a specific area of sample, the system can save its coordinates. If necessary, motorized system can quickly and accurately return to that area, while you are studying another part of the sample

- **Autofocus**

Automated sharpness tuning when working with a sample

- **Remote consultation**

Using the internet you can at any time send images of the specimen to your colleagues for consultation

- **Database**

Digital version of the specimen is stored in the database. You don't need to worry about its safety, loss, physical damage or re-preparation of a sample. You can easily download the digital version of the specimen from the archive

- **Various applications**

- medicine and biology
- cytology and histology
- forensics and pathology
- botany
- osteology
- mycology
- remote consultations and telemedicine
- training for specialists from different fields

WEST MEDICA  
Hegelgasse 19, A-1010,  
Vienna, Austria  
tel.: +43 (1) 804 81 84  
fax: +43 (1) 804 81 85  
office@westmedica.at

WEST MEDICA  
100 King St., Suite 5700,  
Toronto, Ontario, M5X 1C7  
tel.: +1 (416) 915-4245  
fax: +1 (647) 344-6327  
office@westmedica.com

**[www.vision-at.com](http://www.vision-at.com)**  
**[www.westmedica.com](http://www.westmedica.com)**

We reserve the right to change specification without notice.

Distributor

