

# Liquid-based cytology & Vision Cyto<sup>®</sup> Pap

## Workflow in the laboratory

1

### Sample collection



Use specialized instruments and solutions for liquid-based cytology.

2

### Slide preparation



Use automated LBC system.



Use cytological centrifuge.

3

### Slide staining



Slide staining according to PAP protocol.

4

### Slide scanning



Place slides in **Vision Cyto<sup>®</sup> Pap** system for analysis. System automatically scans the slides.



Cells are automatically identified and pre-classified by Bethesda system.

