

∴ **HEINE mini 3000®**

Compact pocket instruments – Maximum quality with minimum dimensions.



➤ **HEINE mini 3000®**

The complete range of high-quality diagnostic instruments in a modern and ergonomic design.



**MINI TO THE MAX.** The handy “mini” series is the name of the compact range of HEINE pocket instruments with the high-quality features of HEINE’s “full-size” instruments. Characteristics include an innovative and compact design, coupled with high-quality, full-range features for the maximum in diagnostic reliability. The series also features automatic power shut-off when returned to the pocket, together with superb stability and manufacturing.

**CONVENIENTLY-SIZED  
AND COMPACT.**

Bigger doesn’t always mean better. In the case of the mini 3000 series, the unique ergonomics and the high-quality features represent the crucial characteristics of superlative quality. The result is exceptional handling in all investigation situations.

**COMPLETE AND PRECISE.**

Our most important goal when developing the mini 3000 series was optimised technical sophistication and enhanced performance. We’ve delivered this through innovative features and improved XHL Xenon Halogen Lens bulb technology for better diagnostic quality and reliability.

**MODERN,  
ERGONOMIC DESIGN.**

Quality meets design – and with the mini 3000 series, the phrase couldn’t be more true. Alongside a modern design, special plastics ensure supreme sturdiness yet low weight. The “mini” series also satisfies every demand in terms of functionality and performance.



HEINE mini 3000®  
BATTERY HANDLE

HEINE mini 3000®  
OPHTHALMOSCOPE

HEINE mini 3000®  
FOCALUX

HEINE mini 3000®  
TONGUE-BLADE HOLDER

Compact battery handle  
suitable for all mini 3000 instruments.



Innovative, compact battery handle in a high-quality, modern design.

- ✦ **High-quality handle: Chrome-finish upper section/refined plastic.** Shockproof, sturdy, non-slip.
- ✦ **Oval, compact, stylish shape.** Lies ergonomically in the hand.
- ✦ **Attachment clip with integrated on/off switch.** Secure. Switches off automatically when returned to the pocket.
- ✦ **Screw-on end cap.** For easy battery change.
- ✦ **Threaded connector made of metal.** Prevents wear.
- ✦ **Complete range of mini 3000 instruments.** Versatile range.
- ✦ **Replaceable batteries.** Size AA.

Ophthalmoscope in a unique, compact design with outstanding optics for even the smallest dimensions. Perfect when on the move and doing home visits.

- ✦ **Enhanced XHL Xenon Halogen Technology.** Concentrated white light for perfect illumination and a brilliant image.
- ✦ **5 different apertures.** A full-function instrument incl. fixation star.
- ✦ **Range of lenses:** -20 to +20D.

Focusing lamp for the examination of the anterior chamber in focal light and for indirect ophthalmoscopy.

- ✦ **Small and light.** Particularly ideal for home visits.
- ✦ **Light spot can be focused on the patient's pupil.** Perfect fundus illumination even with non-dilated pupil.
- ✦ **Enhanced XHL Xenon Halogen Technology for more light.** Very bright, concentrated white light with high light intensity for a bright fundus image.
- ✦ **+3D lens attachment.** Aids accommodation and magnifies image.

Tongue-blade holder with blade ejector for HEINE disposable tongue blades.

- ✦ **Enhanced XHL Xenon Halogen Technology.** Concentrated white light. Concentrated illumination of oral cavity and pharynx.
- ✦ **Finger-tip pressure releases the used blade.** No need to touch the blade after use.
- ✦ **Solid metal construction.** Durable and tough.
- ✦ **Only for HEINE Disposable Tongue-Blades.**

HEINE mini 3000®  
COMBI LAMP



Combination pocket lamp and tongue blade. Durable, multi-purpose instrument.

- ✦ **Many applications.** Can be used as a general illuminator and as an illuminated tongue-blade holder for wooden blades and HEINE disposable blades.
- ✦ **Enhanced XHL Xenon Halogen Lens-lamp with 100% more light.** Bright, focused light.

HEINE mini 3000®  
FIBER OPTIC (F.O.)  
OTOSCOPE



Otoscope with Fiber optic (F.O.) illumination.

- ✦ **Distal Fiber optic (F.O.) illumination with enhanced XHL Xenon Halogen Technology producing 100% more light.** Homogenous, reflex-free illumination of ear canal and tympanum.
- ✦ **Viewing window with 3x magnification** and optimised casing surface for razor-sharp images and minimal reflection. Swivelling viewing window built into instrument.
- ✦ **Swivelling viewing window built into instrument.** Useful for instrumentation, cannot be mislaid.
- ✦ **Insufflation port.** For pneumatic testing of tympanic mobility. Please order connector separately.

HEINE mini 3000®  
OTOSCOPE



Otoscope with direct illumination.

- ✦ **Enhanced XHL Xenon Halogen Technology with 100% more light.** Very bright, concentrated white light for perfect illumination.
- ✦ **Viewing window with 3x magnification** and optimised casing surface for razor-sharp images and minimal reflection.
- ✦ **Swivelling viewing window built into instrument.** Useful for instrumentation, cannot be mislaid.

HEINE mini 3000®  
LARYNGEAL MIRROR



Sturdy, all-purpose laryngeal mirror for optimum illumination of oral cavity and pharynx.

- ✦ **Enhanced XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- ✦ **Solid metal construction.** Durable and tough.

HEINE mini 3000®  
DERMATOSCOPE



HEINE mini 3000®  
CLIP LAMP



HEINE mini-c®  
CLIP LAMP



Dermatoscope with top-quality optics. 10x magnification with focusing optics.

- ✦ **10x magnification.** Large, high-resolution image.
- ✦ **Enhanced XHL Xenon Halogen Technology for more light.** Very bright, concentrated white light for perfect illumination.
- ✦ **Detailed compendium with photographic atlas included.** Simplifies the differential diagnosis of pigmented lesions.

Sturdy, long-life pocket diagnostic lights in a compact, modern design with chrome head. Bright, concentrated light thanks to novel XHL Xenon Halogen lens bulb.

- ✦ **Compact, oval, attractive light.** Lies ergonomically in the hand.
- ✦ **Enhanced XHL Xenon Halogen lens bulb with 100 % more light.** Bright, focused light.

Sturdy, long-life pocket diagnostic light in a compact, modern design with chrome head. mini-c clip lamp only fits mini-c handles.

- ✦ **Compact, oval, attractive light.** Lies ergonomically in the hand.
- ✦ **Enhanced XHL Xenon Halogen lens bulb with 100 % more light.** Bright, focused light.
- ✦ **High-quality handle: Chrome-finish upper section/refined plastic.** Shockproof, sturdy, non-slip.
- ✦ **Switches off automatically when replaced in the pocket.**
- ✦ **Replaceable batteries.** Size AAA.

<b>mini 3000</b>	2.5V
<b>mini 3000 F.O. Otoscope</b> with mini 3000 battery handle and 1 set (4 pcs.) reusable tips and 5 of both 2.5 and 4 mm Ø disposable tips	D-01.70.120
mini 3000 F.O. Otoscope with mini 3000 battery handle and 5 of both 2.5 and 4 mm Ø disposable tips	D-01.70.110
mini 3000 F.O. Otoscope with bulb, without handle and without accessories	D-01.70.106
Insufflation port connector	D-00.80.101
XHL Xenon Halogen spare bulb	X-01.88.105
<b>mini 3000 Otoscope</b> with mini 3000 battery handle, 1 set (4 pcs.) reusable tips and 5 of both 2.5 and 4 mm Ø disposable tips	D-01.70.220
mini 3000 Otoscope with mini 3000 battery handle and 5 of both 2.5 and 4 mm Ø disposable tips	D-01.70.210
mini 3000 Otoscope with bulb, without handle and without accessories	D-01.70.206
XHL Xenon Halogen spare bulb	X-01.88.110
<b>mini 3000 Laryngeal Mirror</b> with plane mirror 20 mm Ø, without handle	D-01.77.101
XHL Xenon Halogen spare bulb	X-01.88.057
<b>mini 3000 Tongue-Blade Holder</b> complete with mini 3000 battery handle and 5 disposable blades	D-01.74.118
Tongue-Blade Holder head with bulb, without handle	D-01.74.100
XHL Xenon Halogen spare bulb	X-01.88.037
<b>mini-c clip lamp</b> in blister pack with 2 AAA batteries and mini-c handle	D-01.73.109
XHL Xenon Halogen spare bulb	X-01.88.108
<b>mini 3000 clip lamp</b> complete with mini 3000 battery handle	D-01.73.131
mini 3000 clip lamp head for connection to mini 3000 battery handle	D-01.73.130
XHL Xenon Halogen spare bulb	X-01.88.107
<b>mini 3000 combi lamp</b> on mini 3000 battery handle with 5 disposable blades	D-01.76.120
mini 3000 combi lamp with bulb, without handle	D-01.76.101
XHL Xenon Halogen spare bulb	X-01.88.107
<b>mini 3000 Ophthalmoscope</b> with mini 3000 battery handle	D-01.71.120
mini 3000 Ophthalmoscope without handle	D-01.71.105
XHL Xenon Halogen spare bulb	X-01.88.106
<b>mini 3000 Focalux</b> with +3D lens with mini 3000 battery handle	D-01.72.141
mini 3000 Focalux without lens attachment with mini 3000 battery handle	D-01.72.139
mini 3000 Focalux without lens, without mini 3000 battery handle	D-01.72.100
Swivel lens attachment +3D only	C-00.14.109
XHL Xenon Halogen spare bulb	X-01.88.042
<b>mini 3000 Dermatoscope</b> with contact plate, mini 3000 battery handle and 10 ml dermatoscopy oil	D-01.78.107
mini 3000 Dermatoscope with contact plate, with bulb, without handle	D-01.78.106
XHL Xenon Halogen spare bulb	X-01.88.109
<b>mini 3000 Battery Handle</b> with two dry batteries IEC LR6 size AA	D-01.79.021
<b>mini 3000 Sets</b>	2.5V
<b>mini 3000 ENT Diagnostic Sets</b>	
Set D-851 complete with mini 3000 Fiber Optic (F.O.) Otoscope, 4 reusable tips, 10 AllSpec disposable tips in hard case with: mini 3000 battery handle with batteries	D-851.10.021
as above, but with mini 3000 Otoscope	D-851.20.021
Set D-861 complete with mini 3000 Fiber Optic (F.O.) Otoscope, 1 set = 4 reusable tips, 10 AllSpec disposable tips in soft pouch with: mini 3000 battery handle with batteries	D-861.10.021
as above, but with mini 3000 Otoscope	D-861.20.021
<b>mini 3000 Ophthalmological Diagnostic Sets</b>	
Set D-852 complete with mini 3000 Ophthalmoscope in hard case, with mini 3000 battery handle with batteries	D-852.10.021
Set D-862 complete with mini 3000 Ophthalmoscope in soft pouch, with mini 3000 battery handle with batteries	D-862.10.021
<b>mini 3000 Combined Diagnostic Sets</b>	
Set D-854 complete with mini 3000 Ophthalmoscope, mini 3000 Otoscope, mini 3000 combi lamp, mini 3000 straight Laryngeal Mirror, straight with plane and concave mirror, 1 set = 4 reusable tips, in hard case, with mini 3000 battery handle with batteries	D-854.10.021
as D-854, but with nose speculum, 10 mm Ø	D-853.10.021
Set D-873 complete with mini 3000 Ophthalmoscope, mini 3000 F.O. Otoscope, 1 set = 4 reusable tips, 5 each of 2.5 and 4 mm Ø AllSpec disposable tips, in hard case, with 2 mini 3000 battery handles with batteries	D-873.11.021
as D-873, but with mini 3000 Otoscope	D-873.21.021
Set D-883 complete with mini 3000 Ophthalmoscope, mini 3000 F.O. Otoscope, 1 set = 4 reusable tips, 5 each of 2.5 and 4 mm Ø AllSpec disposable tips, in soft pouch, with 2 mini 3000 battery handles with batteries	D-883.11.021
as D-883, but with mini 3000 Otoscope	D-883.21.021



## HEINE QUALITY

MADE IN GERMANY

**A company makes history and gives a name to outstanding product quality: HEINE.**

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

We reserve the right to change specifications without notice.

Your distributor

11/06. med 4906 e

## HEINE

HEINE Optotechnik GmbH & Co. KG · Kientalstr. 7 · 82211 Herrsching · Germany  
Tel. +49 (0) 8152/38-0 · Fax +49 (0) 8152/38-202 · E-Mail: [info@heine.com](mailto:info@heine.com) · [www.heine.com](http://www.heine.com)

