

User-friendly, Affordable

Stereo Microscope

EM-30 SERIES

- Robust metal body with rack & pinion
Fresh design, one-point blue parts on white color body
- Light weight, compact body, easy to carry with hand grip
- Stylish, low-profile light base
- Long-life and electricity saving LED light equipped with all models
Transmitted and incident illumination equipped shown below.



• 2 year Warranty

EM-32 Model



EM-30 Series Ordering Information

MODEL	EM-30	EM-31	EM-32	EM-33
Magnification	10X, 30X	20X, 40X	7X - 45X	
Eyepieces	Paired HWF10X/20mm			
Objectives	Turret 1X/3X	Turret 2X/4X	Zoom 0.7X-4.5X	
Viewing head	Binocular 45 degree inclined			Trino. 45 deg.inclined
Diopter	Dioptric adjustment on both eyetubes, Interpupillary 55-75mm			
Stand	Rack & pinion stand with Hand grip			
Focusing block	Coarse adjustment with tension control			
Working distance	100 mm			
Illumination	3W LED incident light/transmitted light, intensity control,			
Accessories	Paired eyeshields, AC electric cord, dust cover			

USA Branch: URL: <http://www.meijitechno.com>
Meiji Techno America
 5895 Rue Ferrari, San Jose, Ca. 95138
 Phone: 408-226-3454
 Fax: 408-226-0900
 Email: info@meijitechno.com



Manufacturer
 322-1, Chikumazawa, Miyoshi-machi,
 Iruma-gun, Saitama 354-0043, Japan
 Phone : +81-49-259-0111
 Fax : +81-49-259-0113
 Email : meiji@meijitechno.co.jp
 URL: <http://www.meijitechno.com>