

HD620-M 11.6 inch Tablet



Specifications:

Model Name	HD620-M
Monitor	11.6 inch HD 1080P LCD monitor
Sensor size	1/3" sensor
Effective pixel	1920 (H) × 1080 (V) Color RETINA monitor
Resolution	60fps@1920*1080, 2.0Mega pixel
Output interface	USB2.0, SD Card, HDMI
Pixel Size	2.80μm × 2.80μm
Power	DC-12V
UI Desktop	Easy mouse control
Mark function	Point coordinates, Reticle, Coordinates, Text explanation
Length measurement	Straight length, broken line length, curve line length, distance between parallel lines, point-line distance
Geometry area	Polygon, Square
Geometry measurement	Line length, radius setting circle, two points setting circle, three points setting circle, concentric circles
Center distance of circles	Radius setting circle, two points setting circle, three points setting circle
Image capture	Support 2M static image
Measurement tool	S-EYE software

More Pictures:

