

## HD750-MTS 9.7 inch Touch Screen Tablet



**Specifications:**

<b>Model Name</b>	HD750-MTS
<b>Category</b>	9.7 inch all in one monitor Touch screen display with Android system 2048 (H) × 1536 (V) resolution
<b>Sensor size</b>	1/2.5 inch
<b>Effective pixel on camera system</b>	5.0Mega pixel
<b>Output interface</b>	USB 2.0, HDMI, TF Card, WIFI
<b>Pixel Size</b>	2.8μm × 2.8μm
<b>Power</b>	DC-12V
<b>UI Desktop</b>	WIFI-mouse operation, fully humanized settings
<b>Indication function</b>	point coordinates, cross line coordinates, text notes
<b>Measurement</b>	Measure line length, curve length, parallel line distance, point line distance
<b>Geometric measurement</b>	radius of the circle line, two round, three point circle, concentric circle
<b>Geometric area</b>	Geometric area polygon and square
<b>Snapshot</b>	Support 5 million static camera snapshot
<b>Video</b>	Video support 1080P video
<b>System</b>	S-eye software

**More Pictures:**

