

HD816-MTS 11.6 inch Touch Screen Tablet



Specifications:

Model Name	Scopepad-LX116W
Category	11.6 inch touch screen LCD monitor
Sensor size	MN34120 1/2.3"
Effective pixel	16Mega pixel resolution, equal as 4032 × 3024 resolution
Resolution	30fps@1080P
Output interface	USB 2.0, SD Card, HDMI, WIFI
Pixel Size	1.335μm × 1.335μm
Power	DC-12V
UI Desktop	Digital design, touch screen or mouse operation, humanized setting
Mark function	Point coordinates, cross lines, coordinate systems, text annotations
Length measuring	Line length, line length, length of curve, parallel line distance and point line distance
Center distance of circles	Radius fixed circle, two point circle, three point fixed circle
Geometry measurement	Line length, radius fixed circle, two point circle, three point fixed circle, concentric circle
Geometry area	Polygon and square
Measuring tool	S-eye software

More Pictures:

