



# MA522

## Super Widefield Eyepiece 10X (Paired) High Eyepoint



MA522 Super Widefield Eyepiece High Eye Point 10X (for EMZ-5H, EMZ-5TRH, EMZ-8TRH, EMZ-10H, EMZ-13H, EMZ-13TRH) Only)

### KEY FEATURES

- 10X High Eyepoint Eyepieces
- Field No. 23.0
- Focal Length: 50.0
- 30.0mm Outer Diameter
- Accepts 25mm reticle
- Made in Japan
- Limited Lifetime Warranty



A pair of optional rubber eyeshields are available for EMZ, EMT and EMF eyepieces (MA906).

### SPECIFICATIONS

#### MA522 (Paired)

A Super Widefield Eyepiece, high eye point 10X (paired) (With 25mm reticle mount) (Fits EMZ-5H, EMZ-5TRH, EMZ-8TRH, EMZ-10H, EMZ-13H, EMZ-13TRH models only)

**NOTE:** Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

#### Optional Eyepieces Available

- **MA501:** Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- **MA718:** Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- **MA502:** Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- **MA520:** Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- **MA503:** Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- **MA535:** Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- **MA504:** Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)
- **MA521:** Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- **MA519:** Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)
- **MA600:** Rubber Eyeshields for SWF eyepieces (paired) (Fits 30.5mm O.D eyepieces) (fits EMZ, EMT, and EMF models)

### PACKAGING INFORMATION

Ships in 8" x 8" x 5" package  
Weight: 0.3 lbs