

## **MX 05**

Economical monocular  
microscope

- Triple objective nosepiece
- 45° inclined head
- 3 objectives achromat: 4x/0.10, 10x/0.25, 40x/0.65
- Widefield eyepieces: 10x/16 mm and 25x
- Built-in LED illumination adjustable 5 V, 1 W
- Up and down illumination



# Specifications

**MICROOPTIX**

## MX 05 monocular microscope

### Standard set

<b>Viewing Head</b>	monocular head, Inclined at 45°
<b>Eyepiece</b>	WF10x/16 and WF25x
<b>Objective</b>	Achromat: 4x/0.10, 10x/0.25, 40x/0.60
<b>Condenser</b>	NA0.65 with Disc Diaphragm
<b>Nosepiece</b>	triple nosepiece
<b>Stage</b>	Plain Stage with slide clips 95x95 mm
<b>Illumination</b>	up and down LED illumination
<b>Focusing System</b>	coaxial coarse and fine adjustment
<b>Additionally (on request)</b>	objectives, digital cameras, software for management of digital albums

## Ordering Information

### Description

### Code

MX 05 monocular microscope, standard set

09.0005.01



We reserve the right to change specification without notice.

West Medica  
Brown-Boveri-Strae 6, B17-1  
2351 Wiener Neudorf, sterreich  
tel.: +43 (0) 2236 892465  
fax: +43 (0) 2236 892464  
vienna@westmedica.com

[www.microoptix.com](http://www.microoptix.com)  
[www.westmedica.com](http://www.westmedica.com)

