

## **MX 1400**

Professional  
stereo microscope

- Binocular ZOOM stereo microscope
- Parallel optical ZOOM system
- Widefield eyepieces 10x/22 mm
- Zoom range: 8–80x
- Plan achromat objectives
- LED illumination adjustable:  
incident and transmitted light



# Specifications

## MX 1400 professional stereo microscope

### Standard set

<b>Optical system</b>	parallel optical ZOOM
<b>Viewing head</b>	binocular Head, 20° Inclination
<b>Interpupillary distance</b>	55–75 mm
<b>Eyepiece</b>	widefield EW 10x/22
<b>ZOOM objective</b>	— 0.8x–8x, ZOOM ratio 1:10 — Plan Achromatic Objective 1x
<b>Working distance</b>	78 mm
<b>Focusing</b>	— coaxial coarse and fine focusing unit — focusing Range 105 mm
<b>Illumination</b>	transmission or Reflection LED Illumination, Brightness Adjustable
<b>Additionally (on request)</b>	eyepieces, lamps, objectives, stands, photo/video trinocular head

## Ordering Information

### Description

### Code

MX 1400 professional stereo microscope, standard set

09.1400.02



We reserve the right to change specification without notice.