

MX 1400

Professional
stereo microscope

- Binocular ZOOM stereo microscope
- Parallel optical ZOOM system
- Widefield eyepieces 10x/22 mm
- Zoom range: 8–80x
- Plan achromat objectives
- LED illumination adjustable:
incident and transmitted light



MX 1400 professional stereo microscope

Standard set

Optical system	parallel optical ZOOM
Viewing head	binocular Head, 20° Inclination
Interpupillary distance	55–75 mm
Eyepiece	widefield EW 10x/22
ZOOM objective	— 0.8–80x, ZOOM ratio 1:10 — Plan Achromatic Objective 1x
Working distance	78 mm
Focusing	— coaxial coarse and fine focusing unit — focusing Range 105 mm
Illumination	transmission or Reflection LED Illumination 5 V, 1 W, Brightness Adjustable
Additionally (on request)	eyepieces, lamps, objectives, stands, photo/video trinocular head

Ordering Information

Description

Code

MX 1400 professional stereo microscope, set 1	09.1400.01
MX 1400 professional stereo microscope, set 2	09.1400.02
MX 1400 professional stereo microscope, set 3	09.1400.03
MX 1400 professional stereo microscope, set 4	09.1400.04
MX 1400 professional stereo microscope, set 5	09.1400.05



We reserve the right to change specification without notice.